



Date: April 10, 2023

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSP 2016

Subject: Microsoft Win Pro Embedded 10 LTSP 2016 Patch release March 2023

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB's validated for March 2023:

KB890830	KB5023697
KB5023788	KB5021045

List of all patches/KB's validated on the Apex System:

KB890830	KB5019964	KB5014702
KB5023788	KB5018411	KB5014164
KB5023697	KB5017396	KB5014026
KB5022838	KB5017305	KB5013952
KB5022289	KB5017095	KB5012596
KB5021235	KB5016884	KB5012170
KB5021045	KB5016622	KB5011570
KB5020627	KB5016058	KB5011495
KB5019966	KB5015808	KB5010359

KB5009557	KB4577586	KB4132216
KB5008873	KB4577015	KB4052908
KB5008207	KB4576750	KB4022619
KB5007192	KB4571694	KB4018073
KB5006669	KB4565912	KB2565063
KB5005698	KB4565511	KB2467173
KB5005573	KB4562561	
KB5005112	KB4561616	
KB5005043	KB4561600	
KB5004238	KB4556813	
KB5003638	KB4550994	
KB5003197	KB4550929	
KB5002323	KB4540723	
KB5002223	KB4540670	
KB5001402	KB4537764	
KB5001347	KB4537759	
KB5001078	KB4535680	
KB5000803	KB4535288	
KB5000685	KB4534271	
KB5000645	KB4530689	
KB4601392	KB4520724	
KB4601318	KB4516115	
KB4598243	KB4516044	
KB4593226	KB4512574	
KB4589210	KB4505422	
KB4589208	KB4503537	
KB4586830	KB4503308	
KB4583465	KB4503267	
KB4583462	KB4498947	
KB4580346	KB4494175	
KB4580325	KB4346087	