



Date: April 10, 2023

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release March 2023

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

List of patches/KB's validated for March 2023:

| | | |
|----------|-----------|-----------|
| KB890830 | KB5023702 | KB5021126 |
|----------|-----------|-----------|

List of all patches/KB's validated on the SecurView System:

| | | |
|-----------|-----------|-----------|
| KB890830 | KB5018419 | KB5011503 |
| KB5023702 | KB5017315 | KB5010351 |
| KB5022840 | KB5016623 | KB5009557 |
| KB5022511 | KB5015811 | KB5008873 |
| KB5022286 | KB5014692 | KB5008218 |
| KB5021237 | KB5013941 | KB5007206 |
| KB5021126 | KB5013641 | KB5005112 |
| KB5020866 | KB5012647 | KB4589208 |
| KB5020627 | KB5012170 | KB4577586 |
| KB5019966 | KB5012128 | KB4535680 |