



**Date:** April 16, 2026

**Author:** Software Engineering

**System:** AWM

**Subsystem:** Windows Server 2019 & Software Security Patches

**Subject:** Windows Server 2019 & Software Security Patch release March 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the AWM System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for March 2026:**

KB5068402

KB5078752

**List of all patches/KB's validated on the AWM System:**

KB5078752	KB5065225	KB5049608
KB5075904	KB5068791	KB5063877
KB5073723	KB5066143	KB5062557
KB5071544	KB5066586	KB5060531

KB5058392	KB5029376	KB5012596
KB5055519	KB5029247	KB5012170
KB5050533	KB5028168	KB5011570
KB890830	KB5027222	KB5011495
KB5053596	KB5027131	KB5010359
KB5052000	KB5026362	KB5009557
KB5050008	KB5025229	KB5009546
KB5049608	KB5023702	KB5008873
KB5048661	KB5022840	KB5008207
KB5046615	KB5022511	KB5007192
KB5044277	KB5022286	KB5006669
KB5044016	KB5021237	KB5005698
KB5043050	KB5021126	KB5005573
KB5042215	KB5020866	KB5005112
KB5041578	KB5020627	KB5005043
KB5040430	KB5019966	KB5004238
KB5039879	KB5019964	KB5003638
KB5039217	KB5018419	KB5003197
KB5037765	KB5018411	KB5001402
KB5036896	KB5017396	KB5001347
KB5036604	KB5017305	KB5001078
KB5035849	KB5017095	KB5000803
KB5034768	KB5016884	KB5000685
KB5034127	KB5016622	KB4601392
KB5033904	KB5016058	KB4601318
KB5033371	KB5015808	KB4598243
KB5032196	KB5014702	KB4593226
KB5031984	KB5014553	KB4589210
KB5031361	KB5014026	KB4589208
KB5030214	KB5013952	KB4586830
KB5029931	KB5013641	KB4580346

KB4577586  
KB4577015  
KB4576750  
KB4571694  
KB4565912  
KB4565511  
KB4562561  
KB4561616  
KB4557397  
KB4551853  
KB4550929  
KB4549947  
KB4540723  
KB4539571  
KB4537759  
KB4535680  
KB4514366  
KB4505225  
KB4494175  
KB4494174  
KB4132216  
KB2565063  
KB2467173