



**Date:** April 14, 2026

**Author:** Software Engineering

**System:** TRIDENT HD

**Subsystem:** Windows 10 & Software Security Patches

**Subject:** Windows 10 & Software Security Patch release March 2026

**Purpose:**

To release the list of validated Microsoft critical patches on the TRIDENT HD System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for March 2026:**

KB5068402

KB5078938

**List of all patches/KB's validated on the TRIDENT HD System:**

KB5068402	KB5073722	KB5066836
KB5078938	KB5071543	KB5065427
KB5075902	KB5065225	KB5065687
KB5075999	KB5068864	KB5063871
KB5073447	KB5070247	KB5062560
	KB5066584	KB5062799

KB5060954	KB5033373	KB5013952
KB5061010	KB5032391	KB5012596
KB5055661	KB5032197	KB5012170
KB5055521	KB5031362	KB5011570
KB5050533	KB5030504	KB5011495
KB890830	KB5030213	KB5010359
KB5054006	KB5029376	KB5009546
KB5053594	KB5029242	KB5008207
KB5052006	KB5028169	KB5007192
KB5050109	KB5027219	KB5006669
KB5049993	KB5026363	KB5005698
KB5048671	KB5025228	KB5005573
KB5046612	KB5023788	KB5005043
KB5044293	KB5023697	KB5004238
KB5043124	KB5022838	KB5003638
KB5043051	KB5022289	KB5003197
KB5042215	KB5021235	KB5002464
KB5041773	KB5021126	KB5002406
KB5041576	KB5021045	KB5001402
KB5040562	KB5019964	KB5001347
KB5040434	KB5018411	KB5001078
KB5039334	KB5017396	KB5000803
KB5039214	KB5017305	KB5000685
KB5037763	KB5017095	KB4601392
KB5037016	KB5016884	KB4601318
KB5036899	KB5016622	KB4598243
KB5035962	KB5016058	KB4593226
KB5035855	KB5015808	KB4589210
KB5034862	KB5014702	KB4586830
KB5034767	KB5014553	KB4583465
KB5034119	KB5014026	KB4580346

KB4577586	KB4537764	KB4498947
KB4577015	KB4535680	KB4494440
KB4576750	KB4534271	KB4494175
KB4571694	KB4530689	KB4493470
KB4565912	KB4525236	KB4489882
KB4565511	KB4521858	KB4487026
KB4562561	KB4520724	KB4485447
KB4561616	KB4519998	KB4346087
KB4557397	KB4516044	KB4132216
KB4556813	KB4512574	KB4091664
KB4550994	KB4512517	KB4022619
KB4550929	KB4509091	KB4018073
KB4540723	KB4507460	KB2565063
KB4540670	KB4505225	KB2467173