



**Date:** 06/02/2020

**Author:** Software Engineering

**System:** Apex

**Subsystem:** Microsoft Windows OS 10

**Subject:** Microsoft Windows 10 OS Patch release May 2020

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Apex systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

<b>System Type</b>	<b>Software Versions</b>
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

**List of patches/KB’s validated for May 2020:**

KB4556813

**List of all patches/KB’s validated on the Apex System:**

KB4556813	KB4505422	KB4018073
KB4550994	KB4494175	KB2565063
KB4540723	KB4132216	KB2467173
KB4537759	KB4022619	