



Date: *May 28, 2021*

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2019 & Software Security Patches

Subject: Windows Server 2019 & Software Security Patch release May 2021

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Note: Install SFW-03452 twice to ensure all patches are applied. The first execution of the patch disk installs the newly required Service Stack and a second run of the patch disk is required to install the latest cumulative patch for May

List of patches/KB's validated for May 2021:

KB5003243 KB5003171 KB5001407

List of all patches/KB's validated on the AWM System:

KB5003243	KB4586793	KB4538461
KB5003171	KB4580325	KB4537759
KB5001407	KB4578966	KB4535680
KB5001404	KB4577586	KB4534273
KB5001342	KB4570720	KB4532691
KB5000859	KB4566424	KB4530715
KB5000822	KB4565349	KB4523205
KB5000685	KB4562562	KB4523204
KB4601393	KB4561608	KB4521862
KB4601345	KB4561600	KB4519338
KB4601060	KB4557397	KB4516115
KB4598480	KB4552924	KB4514366
KB4598230	KB4551853	KB4512578
KB4592440	KB4549949	KB4512577
KB4589208	KB4549947	KB4512517
KB4587735	KB4539571	KB4509091

KB4507460

KB4505225

KB4503537

KB4503267

KB4498947

KB4494440

KB4494174

KB4493470

KB4489882

KB4485447

KB4346087

KB4132216

KB4091664

KB2565063

KB2467173