



Date: May 28, 2021

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release May 2021

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Note: Install SFW-03452 twice to ensure all patches are applied. The first execution of the patch disk installs the newly required Service Stack and a second run of the patch disk is required to install the latest cumulative patch for May

List of patches/KB's validated for May 2021:

KB5003243	KB5003171	KB5001407
-----------	-----------	-----------

List of all patches/KB's validated on the SELENIA DIMENSIONS / 3DIMENSIONS System:

KB890830	KB4601393	KB4580325
KB5003243	KB4601345	KB4577668
KB5003171	KB4601060	KB4577667
KB5001407	KB4598480	KB4577586
KB5001404	KB4598230	KB4570333
KB5001342	KB4592440	KB4570332
KB5000859	KB4589208	KB4566424
KB5000822	KB4587735	KB4565349
KB5000685	KB4586793	KB4562562

KB4561608

KB4551853

KB4514366

KB4561600

KB4549949

KB4505225

KB4558998

KB4549947

KB4494174

KB4558997

KB4539571

KB2565063

KB4557397

KB4537759

KB2467173

KB4552924

KB4535680