



**Date:** June 6, 2022

**Author:** Software Engineering

**System:** TRIDENT

**Subsystem:** Windows 7 & Software Security Patches

**Subject:** Windows 7 & Software Security Patch release May 2022

**Purpose:**

To release the list of validated Microsoft critical patches on the TRIDENT System.

**Note:**

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 004 (Windows 7 ESU License Installation) to the system prior to installation of March 2022 Windows 7 Monthly Patch Disk (SFW-03646) to allow for updates to be installed, due to end of support for Microsoft Windows 7.

**List of patches/KB's validated for May 2022:**

KB5014012	KB5013999	KB5013637
-----------	-----------	-----------

**List of all patches/KB's validated on the TRIDENT System:**

KB976932	KB5012138	KB5009621
KB5014012	KB5012131	KB5008867
KB5013999	KB5011649	KB5008859
KB5013637	KB5011552	KB5008282
KB5012649	KB5011529	KB5007233
KB5012626	KB5010422	KB5006749

KB5006743	KB4592510	KB4561643
KB5006728	KB4592503	KB4556843
KB5005633	KB4592471	KB4556836
KB5005615	KB4586827	KB4556798
KB5005089	KB4586805	KB4555449
KB5005088	KB4585204	KB4552940
KB5004951	KB4580970	KB4552919
KB5004757	KB4580387	KB4550965
KB5004378	KB4580345	KB4541500
KB5004307	KB4578963	KB4538483
KB5004289	KB4578952	KB4537829
KB5004120	KB4577053	KB4537813
KB5003694	KB4577051	KB4536952
KB5003667	KB4576612	KB4534314
KB5003233	KB4575903	KB4534310
KB5003228	KB4571729	KB4533012
KB5001848	KB4571719	KB4532945
KB5001392	KB4570673	KB4532932
KB5001335	KB4569775	KB4531786
KB5000851	KB4569767	KB4530734
KB5000841	KB4565623	KB4530692
KB4601363	KB4565612	KB4525235
KB4601347	KB4565539	KB4525233
KB4600945	KB4565524	KB4523206
KB4598289	KB4565354	KB4523206
KB4598279	KB4562030	KB4520003
KB4597239	KB4561669	KB4519976

KB4516655	KB4483187	KB4038779
KB4516065	KB4480970	KB4034679
KB4516033	KB4480960	KB4025337
KB4515854	KB4480063	KB4022722
KB4512506	KB4480055	KB4022619
KB4512486	KB4474419	KB4020322
KB4511516	KB4471328	KB4019990
KB4507456	KB4468323	KB4019263
KB4507449	KB4467106	KB4015546
KB4507004	KB4462915	KB4012212
KB4506997	KB4457145	KB3212642
KB4505422	KB4457044	KB3205394
KB4503292	KB4457035	KB3197867
KB4503269	KB4343899	KB3192391
KB4499175	KB4339284	KB3177467
KB4499164	KB4338823	KB3172605
KB4495606	KB4284867	KB3161949
KB4495588	KB4103712	KB3159398
KB4493472	KB4096418	KB3156016
KB4493448	KB4093108	KB3150220
KB4490628	KB4074837	KB3139914
KB4489885	KB4074587	KB3139398
KB4489878	KB4057114	KB3138910
KB4486564	KB4056897	KB3138612
KB4486563	KB4054521	KB3133977
KB4483458	KB4048960	KB3127220
KB4483451	KB4041678	KB3126587

KB3124275	KB3035126	KB2894844
KB3122648	KB3030377	KB2893294
KB3110329	KB3023215	KB2888049
KB3109560	KB3022777	KB2884256
KB3109103	KB3021674	KB2882822
KB3108664	KB3019978	KB2871997
KB3108371	KB3018238	KB2868038
KB3101722	KB3011780	KB2864202
KB3097989	KB3010788	KB2862335
KB3092601	KB2992611	KB2862330
KB3086255	KB2991963	KB2862152
KB3080079	KB2984972	KB2861698
KB3078601	KB2979596	KB2840631
KB3075220	KB2978742	KB2834140
KB3074543	KB2977292	KB2813430
KB3072305	KB2973201	KB2813347
KB3071756	KB2973112	KB2807986
KB3067903	KB2972211	KB2789645
KB3061518	KB2972100	KB2786081
KB3060716	KB2968294	KB2770660
KB3059317	KB2943357	KB2758857
KB3055642	KB2937610	KB2742599
KB3046269	KB2931356	KB2736422
KB3046017	KB2912390	KB2731771
KB3045685	KB2911501	KB2729452
KB3042058	KB2900986	KB2729094
KB3037574	KB2894854	KB2727528

KB2719033	KB2506014
KB2716513	KB2491683
KB2705219	KB2479943
KB2698365	KB2467173
KB2690533	
KB2685939	
KB2670838	
KB2667402	
KB2656356	
KB2654428	
KB2653956	
KB2639308	
KB2631813	
KB2621440	
KB2620704	
KB2604115	
KB2585542	
KB2579686	
KB2565063	
KB2564958	
KB2560656	
KB2552343	
KB2533623	
KB2533552	
KB2532531	
KB2510531	
KB2506212	