



Date: June 13, 2023

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSB 2016

Subject: Microsoft Win Pro Embedded 10 LTSB 2016 Patch release May 2023

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSB 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB’s validated for May 2023:

KB5026363

KB5002369

List of all patches/KB’s validated on the Apex System:

KB890830

KB5019966

KB5014702

KB5026363

KB5019964

KB5014164

KB5025228

KB5018411

KB5014026

KB5023788

KB5017396

KB5013952

KB5023697

KB5017305

KB5012596

KB5022838

KB5017095

KB5012170

KB5022289

KB5016884

KB5011570

KB5021235

KB5016622

KB5011495

KB5021045

KB5016058

KB5010359

KB5020627

KB5015808

KB5009557

KB5008873	KB4593226	KB4537759
KB5008207	KB4589210	KB4535680
KB5007192	KB4589208	KB4535288
KB5006669	KB4586830	KB4534271
KB5005698	KB4583465	KB4530689
KB5005573	KB4583462	KB4520724
KB5005112	KB4580346	KB4516115
KB5005043	KB4580325	KB4516044
KB5004238	KB4577586	KB4512574
KB5003638	KB4577015	KB4505422
KB5003197	KB4576750	KB4503537
KB5002369	KB4571694	KB4503308
KB5002323	KB4565912	KB4503267
KB5002223	KB4565511	KB4498947
KB5001402	KB4562561	KB4494175
KB5001347	KB4561616	KB4346087
KB5001078	KB4561600	KB4132216
KB5000803	KB4556813	KB4052908
KB5000685	KB4550994	KB4022619
KB5000645	KB4550929	KB4018073
KB4601392	KB4540723	KB2565063
KB4601318	KB4540670	KB2467173
KB4598243	KB4537764	