

Date: June 09, 2025

Author: Software Engineering

System: InSight

Subsystem: Microsoft Windows 10 Enterprise 2016 LTSC

Subject: Microsoft Windows 10 Enterprise 2016 LTSC Patch release May 2025

Purpose:

To release the list of validated Microsoft critical patches on the InSight systems.

System Type	Software Versions
Workstation	InSight v6.1.4

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for May 2025:

KB5058524

KB5058383

List of all patches/KB's validated on the InSight System:

KB5058524

KB5039885

KB5034862

KB5058383

KB5039334

KB5034767

KB5055661

KB5039214

KB5034119

KB5055521

KB5037763

KB5033910

KB5049614

KB5037016

KB5033904

KB5002702

KB5036899

KB5033373

KB5041773

KB5036609

KB5032391

KB5041576

KB5035962

KB5032197

KB5040562

KB5035855

KB5031989

KB5040434

KB5035849

KB5031362

KB5030504	KB5011570	KB4598243
KB5030213	KB5011495	KB4593226
KB5029924	KB5010359	KB4589210
KB5029242	KB5009546	KB4589208
KB5028169	KB5008207	KB4586830
KB5027219	KB5007192	KB4580346
KB5027123	KB5006669	KB4580325
KB5026363	KB5005698	KB4577586
KB5025228	KB5005573	KB4577015
KB5023788	KB5005112	KB4576750
KB5023697	KB5005043	KB4571694
KB5022838	KB5004424	KB4565912
KB5022503	KB5004238	KB4565511
KB5022289	KB5003638	KB4562561
KB5021235	KB5003197	KB4561616
KB5020873	KB5002662	KB4561600
KB5020614	KB5002619	KB4556813
KB5019964	KB5002542	KB4550994
KB5018411	KB5002520	KB4550929
KB5017396	KB5002497	KB4540723
KB5017305	KB5002464	KB4540670
KB5017095	KB5002406	KB4537764
KB5016622	KB5002369	KB4537759
KB5016058	KB5002323	KB4535680
KB5015808	KB5002223	KB4534271
KB5014702	KB5001402	KB4530689
KB5014026	KB5001347	KB4525236
KB5013952	KB5001078	KB4520724
KB5013625	KB5000803	KB4516115
KB5012596	KB4601392	KB4494175
KB5012170	KB4601318	KB4132216

KB2467173