



**Date:** June 09, 2025

**Author:** Software Engineering

**System:** SecurView Manager – Windows Server 2016

**Subsystem:** Windows 10 Server 2016 & Software Security Patches

**Subject:** Windows 10 Server 2016 & Software Security Patch release May 2025

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2016 System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for May 2025:**

KB5058524                      KB5058383

**List of all patches/KB's validated on the SecurView Manager – Windows Server 2016 System:**

|           |           |           |
|-----------|-----------|-----------|
| KB5058524 | KB5052006 | KB5043124 |
| KB5055661 | KB5050109 | KB5043051 |
| KB5055521 | KB5049993 | KB5041773 |
| KB890830  | KB5048671 | KB5041576 |
| KB5054006 | KB5046612 | KB5040562 |
| KB5053594 | KB5044293 | KB5040434 |

|           |           |           |
|-----------|-----------|-----------|
| KB5039334 | KB5032391 | KB5023697 |
| KB5039214 | KB5032197 | KB5022838 |
| KB5037763 | KB5031362 | KB5022289 |
| KB5037016 | KB5030504 | KB5017396 |
| KB5036899 | KB5030213 | KB5012170 |
| KB5035962 | KB5029242 | KB5005698 |
| KB5035855 | KB5028169 | KB4589210 |
| KB5034862 | KB5027219 | KB4494175 |
| KB5034767 | KB5026363 | KB4132216 |
| KB5034119 | KB5025228 |           |
| KB5033373 | KB5023788 |           |