



Date: June 08, 2026

Author: Software Engineering

System: BREVERA / CORLUMINA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release May 2026

Purpose:

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for May 2026:

KB5087538

KB5087061

List of all patches/KB's validated on the BREVERA / CORLUMINA System:

| | | |
|-----------|-----------|-----------|
| KB5087061 | KB5050008 | KB5031361 |
| KB5087538 | KB5049608 | KB5030214 |
| KB5082123 | KB5048661 | KB5029931 |
| KB5082413 | KB5046615 | KB5029376 |
| KB5068402 | KB5044277 | KB5029247 |
| KB5078752 | KB5044016 | KB5028168 |
| KB5075904 | KB5043050 | KB5027222 |
| KB5073723 | KB5042215 | KB5027131 |
| KB5071544 | KB5041578 | KB5026363 |
| KB5068791 | KB5040430 | KB5026362 |
| KB5065225 | KB5039879 | KB5025229 |
| KB5066143 | KB5039217 | KB5025228 |
| KB5066586 | KB5037765 | KB5023788 |
| KB5065428 | KB5036896 | KB5023697 |
| KB5063877 | KB5036604 | KB5022838 |
| KB5062557 | KB5035849 | KB5022511 |
| KB5060531 | KB5034768 | KB5022289 |
| KB5058392 | KB5034127 | KB5021235 |
| KB5055519 | KB5033904 | KB5021126 |
| KB5050533 | KB5033371 | KB5019964 |
| KB890830 | KB5032196 | KB5018411 |
| KB5053596 | KB5031984 | KB5017396 |
| KB5052000 | | |

| | | |
|-----------|-----------|-----------|
| KB5017305 | KB5001347 | KB4535680 |
| KB5017095 | KB5001078 | KB4534271 |
| KB5016884 | KB5000803 | KB4530689 |
| KB5016622 | KB5000685 | KB4525236 |
| KB5016058 | KB4601392 | KB4521858 |
| KB5015808 | KB4601318 | KB4520724 |
| KB5014702 | KB4598243 | KB4519998 |
| KB5014553 | KB4593226 | KB4516044 |
| KB5014026 | KB4589210 | KB4512574 |
| KB5013952 | KB4589208 | KB4512517 |
| KB5012596 | KB4586830 | KB4509091 |
| KB5012170 | KB4580346 | KB4507460 |
| KB5011570 | KB4577586 | KB4505225 |
| KB5011495 | KB4577015 | KB4498947 |
| KB5010359 | KB4576750 | KB4494440 |
| KB5009546 | KB4571694 | KB4494175 |
| KB5008207 | KB4565912 | KB4494174 |
| KB5007192 | KB4565511 | KB4493470 |
| KB5006669 | KB4562561 | KB4489882 |
| KB5005698 | KB4561616 | KB4487026 |
| KB5005573 | KB4557397 | KB4485447 |
| KB5005112 | KB4556813 | KB4346087 |
| KB5005043 | KB4550994 | KB4132216 |
| KB5004238 | KB4550929 | KB4091664 |
| KB5003638 | KB4540723 | KB2565063 |
| KB5003197 | KB4540670 | KB2467173 |
| KB5001402 | KB4537764 | |

