



**Date:** June 08, 2026

**Author:** Software Engineering

**System:** Cenova

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release May 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Cenova System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for May 2026:**

KB5087538

KB5087061

**List of all patches/KB's validated on the Cenova System:**

KB5087061	KB5050008	KB5030214
KB5087538	KB5049608	KB5029931
KB5082123	KB5048661	KB5029376
KB5082413	KB5046615	KB5029247
KB5068402	KB5044277	KB5028168
KB5078752	KB5044016	KB5027222
KB5075904	KB5043050	KB5027131
KB5073723	KB5042215	KB5027124
KB5071544	KB5041578	KB5026362
KB5065225	KB5040430	KB5025229
KB5068791	KB5039879	KB5023702
KB5066143	KB5039217	KB5022840
KB5066586	KB5037765	KB5022511
KB5065428	KB5036896	KB5022286
KB5063877	KB5036604	KB5021237
KB5062557	KB5035849	KB5021126
KB5060531	KB5034768	KB5020866
KB5058392	KB5034127	KB5020627
KB5055519	KB5033904	KB5019966
KB5050533	KB5033371	KB5018419
KB890830	KB5032196	KB5017315
KB5053596	KB5031984	KB5016884
KB5052000	KB5031361	KB5016623

KB5015811	KB5005568	KB5000685
KB5014692	KB5005112	KB4601393
KB5014553	KB5005030	KB4601345
KB5013941	KB5004244	KB4601060
KB5013641	KB5003711	KB4589210
KB5012647	KB5003646	KB4589208
KB5012170	KB5003638	KB4580325
KB5012128	KB5003243	KB4577586
KB5011503	KB5003171	KB4557397
KB5010351	KB5001407	KB4539571
KB5009557	KB5001404	KB4535680
KB5008873	KB5001402	KB4494174
KB5008218	KB5001342	KB4465065
KB5007206	KB5000859	
KB5006672	KB5000822	