



Date: June 08, 2026

Author: Software Engineering

System: Faxitron Path Vision

Subsystem: Windows 10 LTSC 2021 & Software Security Patches

Subject: Windows 10 LTSC 2021 & Software Security Patch release May 2026

Purpose:

To release the list of validated Microsoft critical patches on the Faxitron Path Vision System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for May 2026:

KB5087544

KB5087064

List of all patches/KB's validated on the Faxitron Path Vision System:

KB5087064

KB5087544

KB5012170

KB5082200

KB5082426

KB5066135

KB5078885