



Date: June 08, 2026

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release May 2026

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

KB5019964	KB5004238	KB4550929
KB5018411	KB5003638	KB4540723
KB5017396	KB5003197	KB4540670
KB5017305	KB5001402	KB4537764
KB5017095	KB5001347	KB4535680
KB5016884	KB5001078	KB4534271
KB5016622	KB5000803	KB4530689
KB5016058	KB5000685	KB4525236
KB5015808	KB4601392	KB4521858
KB5014702	KB4601318	KB4520724
KB5014553	KB4601051	KB4519998
KB5014026	KB4598243	KB4516044
KB5013952	KB4593226	KB4512574
KB5013625	KB4589210	KB4512517
KB5012596	KB4586830	KB4509091
KB5012170	KB4580346	KB4507460
KB5012118	KB4578969	KB4505225
KB5011570	KB4577586	KB4503537
KB5011495	KB4577015	KB4503267
KB5010359	KB4576750	KB4498947
KB5009546	KB4571694	KB4494440
KB5008877	KB4565912	KB4494175
KB5008207	KB4565511	KB4493472
KB5007192	KB4562561	KB4493470
KB5006669	KB4561616	KB4493448
KB5005698	KB4557397	KB4493435
KB5005573	KB4556813	KB4489882
KB5005043	KB4550994	KB4485447

KB4346087

KB4091664

KB2467173

KB4132216

KB2565063