



**Date:** June 08, 2026

**Author:** Software Engineering

**System:** SecurView Manager – Windows Server 2022

**Subsystem:** Windows 10 Server 2022 & Software Security Patches

**Subject:** Windows 10 Server 2022 & Software Security Patch release May 2026

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2022 System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for May 2026:**

KB5087545

KB5087068

**List of all patches/KB's validated on the SecurView Manager – Windows Server 2022 System:**

KB5087068

KB5087545

KB5082427

KB5082142

KB5066139

KB5078766