



**Date:** *January 11, 2021*

**Author:** Software Engineering

**System:** TRIDENT

**Subsystem:** Windows 7 & Software Security Patches

**Subject:** Windows 7 & Software Security Patch release November 2020

---

**Purpose:**

To release the list of validated Microsoft critical patches on the TRIDENT System.

**Note:**

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 001 (Windows 7 ESU License Installation) to the system prior to installation of November 2020 Windows 7 Monthly Patch Disk (SFW-03339) to allow for updates to be installed, due to end of support for Microsoft Windows 7.
- After installation of SFW-03339, please check if KB4586827 is installed. If not installed, please download and install KB4586827 manually.

**List of patches/KB's validated for November 2020:**

KB4586827                      KB4585204

**List of all patches/KB's validated on the TRIDENT System:**

KB976932	KB4580387	KB4577053
KB4586827	KB4580345	KB4577051
KB4585204	KB4578963	KB4576612
KB4580970	KB4578952	KB4575903

KB4571729	KB4533012	KB4499164
KB4571719	KB4532945	KB4495606
KB4570673	KB4532932	KB4495588
KB4569775	KB4531786	KB4493472
KB4569767	KB4530734	KB4493448
KB4565623	KB4530692	KB4490628
KB4565612	KB4525235	KB4489885
KB4565539	KB4525233	KB4489878
KB4565524	KB4523206	KB4486564
KB4565354	KB4523206	KB4486563
KB4562030	KB4520003	KB4483458
KB4561669	KB4519976	KB4483451
KB4561643	KB4516655	KB4483187
KB4556843	KB4516065	KB4480970
KB4556836	KB4516033	KB4480960
KB4556798	KB4515854	KB4480063
KB4555449	KB4512506	KB4480055
KB4552940	KB4512486	KB4474419
KB4552919	KB4511516	KB4471328
KB4550965	KB4507456	KB4468323
KB4541500	KB4507449	KB4467106
KB4538483	KB4507004	KB4462915
KB4537829	KB4506997	KB4457145
KB4537813	KB4505422	KB4457044
KB4536952	KB4503292	KB4457035
KB4534314	KB4503269	KB4343899
KB4534310	KB4499175	KB4339284

KB4338823	KB3172605	KB3072305
KB4284867	KB3161949	KB3071756
KB4103712	KB3159398	KB3067903
KB4096418	KB3156016	KB3061518
KB4093108	KB3150220	KB3060716
KB4074837	KB3139914	KB3059317
KB4074587	KB3139398	KB3055642
KB4057114	KB3138910	KB3046269
KB4056897	KB3138612	KB3046017
KB4054521	KB3133977	KB3045685
KB4048960	KB3127220	KB3042058
KB4041678	KB3126587	KB3037574
KB4038779	KB3124275	KB3035126
KB4034679	KB3122648	KB3030377
KB4025337	KB3110329	KB3023215
KB4022722	KB3109560	KB3022777
KB4022619	KB3109103	KB3021674
KB4020322	KB3108664	KB3019978
KB4019990	KB3108371	KB3018238
KB4019263	KB3101722	KB3011780
KB4015546	KB3097989	KB3010788
KB4012212	KB3092601	KB2992611
KB3212642	KB3086255	KB2991963
KB3205394	KB3080079	KB2984972
KB3197867	KB3078601	KB2979596
KB3192391	KB3075220	KB2978742
KB3177467	KB3074543	KB2977292

KB2973201	KB2813347	KB2604115
KB2973112	KB2807986	KB2585542
KB2972211	KB2789645	KB2579686
KB2972100	KB2786081	KB2565063
KB2968294	KB2770660	KB2564958
KB2943357	KB2758857	KB2560656
KB2937610	KB2742599	KB2552343
KB2931356	KB2736422	KB2533623
KB2912390	KB2731771	KB2533552
KB2911501	KB2729452	KB2532531
KB2900986	KB2729094	KB2510531
KB2894854	KB2727528	KB2506212
KB2894844	KB2719033	KB2506014
KB2893294	KB2716513	KB2491683
KB2888049	KB2705219	KB2479943
KB2884256	KB2698365	KB2467173
KB2882822	KB2690533	
KB2871997	KB2685939	
KB2868038	KB2670838	
KB2864202	KB2667402	
KB2862335	KB2656356	
KB2862330	KB2654428	
KB2862152	KB2653956	
KB2861698	KB2639308	
KB2840631	KB2631813	
KB2834140	KB2621440	
KB2813430	KB2620704	