



Date: *December 7, 2020*

Author: Software Engineering

System: UNIFI WORKSPACE SOFTWARE

Subsystem: Windows 10 LTSC & Software Security Patches

Subject: Windows 10 LTSC & Software Security Patch release November 2020

Purpose:

To release the list of validated Microsoft critical patches on the Unifi Workspace Software System.

List of patches/KB's validated for November 2020:

KB4586830

List of all patches/KB's validated on the Unifi Workspace Software System:

KB4586830	KB4561616	KB4494175
KB4580346	KB4561600	KB4346087
KB4580325	KB4556813	KB4132216
KB4577015	KB4550994	KB4052908
KB4576750	KB4550929	KB2565063
KB4565912	KB4549825	KB2467173
KB4565511	KB4540723	
KB4562561	KB4537759	