



Date: 12/20/2021

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSP 2016

Subject: Microsoft Win Pro Embedded 10 LTSP 2016 Patch release November 2021

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB’s validated for November 2021:

KB5007192 KB4022619
KB4583462 KB4018073

List of all patches/KB’s validated on the Apex System:

KB890830 KB5003197 KB4598243
KB5007192 KB5001402 KB4593226
KB5006669 KB5001347 KB4589210
KB5005698 KB5001078 KB4586830
KB5005573 KB5000803 KB4583462
KB5005043 KB5000645 KB4580346
KB5004238 KB4601392 KB4580325
KB5003638 KB4601318 KB4577586

KB4577015	KB4052908
KB4576750	KB4022619
KB4571694	KB4018073
KB4565912	KB2565063
KB4565511	KB2467173
KB4562561	
KB4561616	
KB4561600	
KB4556813	
KB4550994	
KB4550929	
KB4540723	
KB4540670	
KB4537764	
KB4537759	
KB4535680	
KB4535288	
KB4534271	
KB4530689	
KB4520724	
KB4516115	
KB4516044	
KB4512574	
KB4505422	
KB4503537	
KB4503308	
KB4503267	
KB4498947	
KB4494175	
KB4346087	
KB4132216	