



Date: December 9, 2022

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSP 2016

Subject: Microsoft Win Pro Embedded 10 LTSP 2016 Patch release November 2022

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note 1: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB's validated for November 2022:

KB5020627	KB5005112	KB4589208
KB5019966	KB5002223	

List of all patches/KB's validated on the Apex System:

KB890830	KB5016058	KB5011570
KB5020627	KB5015808	KB5011495
KB5019966	KB5014702	KB5010359
KB5018411	KB5014164	KB5009557
KB5017396	KB5014026	KB5008873
KB5017305	KB5013952	KB5008207
KB5017095	KB5012596	KB5007192
KB5016622	KB5012170	KB5006669

KB5005698	KB4561616
KB5005573	KB4561600
KB5005112	KB4556813
KB5005043	KB4550994
KB5004238	KB4550929
KB5003638	KB4540723
KB5003197	KB4540670
KB5002223	KB4537764
KB5001402	KB4537759
KB5001347	KB4535680
KB5001078	KB4535288
KB5000803	KB4534271
KB5000685	KB4530689
KB5000645	KB4520724
KB4601392	KB4516115
KB4601318	KB4516044
KB4598243	KB4512574
KB4593226	KB4505422
KB4589210	KB4503537
KB4589208	KB4503308
KB4586830	KB4503267
KB4583462	KB4498947
KB4580346	KB4494175
KB4580325	KB4346087
KB4577586	KB4132216
KB4577015	KB4052908
KB4576750	KB4022619
KB4571694	KB4018073
KB4565912	KB2565063
KB4565511	KB2467173
KB4562561	