

Date: December 9, 2022

Author: Software Engineering

System: Apex

Subsystem: Win Pro Embedded 10 LTSC 2019

Subject: Microsoft Win Pro Embedded 10 LTSC 2019 Patch release November 2022

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSC 2019 systems.

List of patches/KB's validated for November 2022:

KB5019964	KB5005698	KB4589210
KB5017396	KB5002223	KB4132216

List of all patches/KB's validated on the Apex System:

KB5019964	KB5012128	KB5003711
KB5018419	KB5011503	KB5002223
KB5017396	KB5010351	KB4589210
KB5017315	KB5009557	KB4589208
KB5016623	KB5008873	KB4583462
KB5015811	KB5008218	KB4577586
KB5014692	KB5007206	KB4535680
KB5014164	KB5006672	KB4494174
KB5013941	KB5005698	KB4132216
KB5013641	KB5005568	KB4022619
KB5012647	KB5005112	KB4018073
KB5012170	KB5005030	