



**Date:** December 19, 2023

**Author:** Software Engineering

**System:** Path Vision

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release November 2023

---

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Path Vision System.

Note: Core Vision is on similar product line as Path Vision and both use the same executable so patches validated on Path Vision system can be installed on Core Visions system as well.

**List of patches/KB's validated for November 2023:**

KB5032196                      KB5031984

**List of all patches/KB's validated on the Path Vision System:**

KB890830	KB5025229	KB5015811
KB5032196	KB5023702	KB5014692
KB5031984	KB5022840	KB5013941
KB5031361	KB5022511	KB5013641
KB5030214	KB5021237	KB5012170
KB5029931	KB5020866	KB5005112
KB5029247	KB5020627	KB4589208
KB5028168	KB5019966	KB4577586
KB5027222	KB5018419	KB4535680
KB5027131	KB5017315	KB4494174
KB5026362	KB5016623	