

Date: December 08, 2025

Author: Software Engineering

System: Apex

Subsystem: Win Pro Embedded 10 LTSC 2019

Subject: Microsoft Win Pro Embedded 10 LTSC 2019 Patch release November 2025

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSC 2019 systems.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2019:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for November 2025:

KB5068791

KB5065222

List of all patches/KB's validated on the Apex System:

KB5065222

KB5062557

KB5053596

KB5068791

KB5060531

KB5052000

KB5066143

KB5058392

KB5050008

KB5066586

KB5055519

KB5049608

KB5065428

KB5054833

KB5048661

KB5063877

KB890830

KB5046615

KB5044277	KB5026362	KB5008218
KB5044016	KB5025229	KB5007206
KB5043050	KB5023702	KB5006672
KB5042749	KB5022840	KB5005698
KB5041578	KB5022511	KB5005568
KB5040430	KB5022286	KB5005112
KB5039879	KB5021237	KB5005030
KB5039217	KB5021045	KB5003711
KB5037765	KB5020866	KB5002542
KB5036896	KB5019964	KB5002520
KB5036604	KB5018419	KB5002497
KB5035849	KB5017396	KB5002464
KB5034768	KB5017315	KB5002406
KB5034127	KB5016623	KB5002369
KB5033904	KB5015811	KB5002323
KB5033688	KB5014692	KB5002223
KB5033371	KB5014164	KB4589210
KB5032196	KB5013941	KB4589208
KB5031984	KB5013641	KB4583465
KB5031361	KB5012647	KB4583462
KB5030214	KB5012170	KB4577586
KB5029931	KB5012128	KB4535680
KB5029247	KB5011503	KB4494174
KB5028168	KB5010351	KB4132216
KB5027222	KB5009557	KB4022619
KB5027131	KB5008873	KB4018073