



Date: December 09, 2025

Author: Software Engineering

System: BREVERA / CORLUMINA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
November

Purpose:

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for November 2025:

KB5068791

KB5065225

List of all patches/KB's validated on the BREVERA / CORLUMINA System:

KB5068791	KB5039217	KB5022838
KB5065225	KB5037765	KB5022511
KB5066143	KB5036896	KB5022289
KB5066586	KB5036604	KB5021235
KB5065428	KB5035849	KB5021126
KB5063877	KB5034768	KB5019964
KB5062557	KB5034127	KB5018411
KB5060531	KB5033904	KB5017396
KB5058392	KB5033371	KB5017305
KB5055519	KB5032196	KB5017095
KB5050533	KB5031984	KB5016884
KB890830	KB5031361	KB5016622
KB5053596	KB5030214	KB5016058
KB5052000	KB5029931	KB5015808
KB5050008	KB5029376	KB5014702
KB5049608	KB5029247	KB5014553
KB5048661	KB5028168	KB5014026
KB5046615	KB5027222	KB5013952
KB5044277	KB5027131	KB5012596
KB5044016	KB5026363	KB5012170
KB5043050	KB5026362	KB5011570
KB5042215	KB5025229	KB5011495
KB5041578	KB5025228	KB5010359
KB5040430	KB5023788	KB5009546
KB5039879	KB5023697	KB5008207

KB5007192	KB4580346	KB4520724
KB5006669	KB4577586	KB4519998
KB5005698	KB4577015	KB4516044
KB5005573	KB4576750	KB4512574
KB5005112	KB4571694	KB4512517
KB5005043	KB4565912	KB4509091
KB5004238	KB4565511	KB4507460
KB5003638	KB4562561	KB4505225
KB5003197	KB4561616	KB4498947
KB5001402	KB4557397	KB4494440
KB5001347	KB4556813	KB4494175
KB5001078	KB4550994	KB4494174
KB5000803	KB4550929	KB4493470
KB5000685	KB4540723	KB4489882
KB4601392	KB4540670	KB4487026
KB4601318	KB4537764	KB4485447
KB4598243	KB4535680	KB4346087
KB4593226	KB4534271	KB4132216
KB4589210	KB4530689	KB4091664
KB4589208	KB4525236	KB2565063
KB4586830	KB4521858	KB2467173