



Date: December 09, 2025

Author: Software Engineering

System: Cenova

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
November 2025

Purpose:

To release the list of validated Microsoft critical patches on the Cenova System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for November 2025:

KB5068791

KB5065225

List of all patches/KB's validated on the Cenova System:

KB5065225	KB5037765	KB5021126
KB5068791	KB5036896	KB5020866
KB5066143	KB5036604	KB5020627
KB5066586	KB5035849	KB5019966
KB5065428	KB5034768	KB5018419
KB5063877	KB5034127	KB5017315
KB5062557	KB5033904	KB5016884
KB5060531	KB5033371	KB5016623
KB5058392	KB5032196	KB5015811
KB5055519	KB5031984	KB5014692
KB5050533	KB5031361	KB5014553
KB890830	KB5030214	KB5013941
KB5053596	KB5029931	KB5013641
KB5052000	KB5029376	KB5012647
KB5050008	KB5029247	KB5012170
KB5049608	KB5028168	KB5012128
KB5048661	KB5027222	KB5011503
KB5046615	KB5027131	KB5010351
KB5044277	KB5027124	KB5009557
KB5044016	KB5026362	KB5008873
KB5043050	KB5025229	KB5008218
KB5042215	KB5023702	KB5007206
KB5041578	KB5022840	KB5006672
KB5040430	KB5022511	KB5005568
KB5039879	KB5022286	KB5005112
KB5039217	KB5021237	KB5005030

KB5004244	KB5001342	KB4580325
KB5003711	KB5000859	KB4577586
KB5003646	KB5000822	KB4557397
KB5003638	KB5000685	KB4539571
KB5003243	KB4601393	KB4535680
KB5003171	KB4601345	KB4494174
KB5001407	KB4601060	KB4465065
KB5001404	KB4589210	
KB5001402	KB4589208	