

Date: December 09, 2025

Author: Software Engineering

System: InSight

Subsystem: Microsoft Windows 10 Enterprise 2019 LTSC

Subject: Microsoft Windows 10 Enterprise 2019 LTSC Patch release November 2025

Purpose:

To release the list of validated Microsoft critical patches on the Insights systems.

System Type	Software Versions
Workstation	InSight v6.1.4

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for November 2025:

KB5068791

List of all patches/KB's validated on the Insight System:

KB5068791	KB5053596	KB5040430
KB5066143	KB5052000	KB5039879
KB5066586	KB5050008	KB5039217
KB5065428	KB5049608	KB5037765
KB5063877	KB5048661	KB5036896
KB5062557	KB5046615	KB5036604
KB5060531	KB5044277	KB5035849
KB5058392	KB5044016	KB5034768
KB5055519	KB5043050	KB5034127
KB890830	KB5041578	KB5033904

KB5033371
KB5032196
KB5031984
KB5031361
KB5030214
KB5029931
KB5029247
KB5028168
KB5027222
KB5027131
KB5026362
KB5025229
KB5023702
KB5022840
KB5022511
KB5012170
KB5005112
KB5004424
KB4589208
KB4577586