



Date: December 09, 2025

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release
November 2025

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for November 2025:

KB5065225

KB5070247

KB5068864

**List of all patches/KB's validated on the SELENIA DIMENSIONS /
3DIMENSIONS System:**

KB5065225	KB5070247	KB5066836
KB5068864	KB5066584	KB5065427
KB5065687	KB5041576	KB5025228
KB5063871	KB5040562	KB5023788
KB5062560	KB5040434	KB5023697
KB5062799	KB5039334	KB5022838
KB5060954	KB5039214	KB5022289
KB5061010	KB5037763	KB5021235
KB5058524	KB5037016	KB5021126
KB5058383	KB5036899	KB5019964
KB5055661	KB5035962	KB5018411
KB5055521	KB5035855	KB5017396
KB5050533	KB5034862	KB5017305
KB890830	KB5034767	KB5017095
KB5054006	KB5034119	KB5016884
KB5053594	KB5033373	KB5016622
KB5052006	KB5032391	KB5016058
KB5050109	KB5032197	KB5015808
KB5049993	KB5031362	KB5014702
KB5048671	KB5030504	KB5014553
KB5046612	KB5030213	KB5014026
KB5044293	KB5029376	KB5013952
KB5043124	KB5029242	KB5013625
KB5043051	KB5028169	KB5012596
KB5042215	KB5027219	KB5012170
KB5041773	KB5026363	KB5012118

KB5011570	KB4578969	KB4507460
KB5011495	KB4577586	KB4505225
KB5010359	KB4577015	KB4503537
KB5009546	KB4576750	KB4503267
KB5008877	KB4571694	KB4498947
KB5008207	KB4565912	KB4494440
KB5007192	KB4565511	KB4494175
KB5006669	KB4562561	KB4493472
KB5005698	KB4561616	KB4493470
KB5005573	KB4557397	KB4493448
KB5005043	KB4556813	KB4493435
KB5004238	KB4550994	KB4489882
KB5003638	KB4550929	KB4485447
KB5003197	KB4540723	KB4346087
KB5001402	KB4540670	KB4132216
KB5001347	KB4537764	KB4091664
KB5001078	KB4535680	KB2565063
KB5000803	KB4534271	KB2467173
KB5000685	KB4530689	
KB4601392	KB4525236	
KB4601318	KB4521858	
KB4601051	KB4520724	
KB4598243	KB4519998	
KB4593226	KB4516044	
KB4589210	KB4512574	
KB4586830	KB4512517	
KB4580346	KB4509091	