



Date: December 09, 2025

Author: Software Engineering

System: SecurView Manager – Windows Server 2016

Subsystem: Windows 10 Server 2016 & Software Security Patches

Subject: Windows 10 Server 2016 & Software Security Patch release November 2025

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2016 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.21 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for November 2025:

KB5070247 KB5068864

List of all patches/KB's validated on the SecurView Manager – Windows Server 2016 System:

KB5068864	KB5066584	KB5065427
KB5070247	KB5066836	KB5065687
KB5063871	KB5061010	KB5055661
KB5062560	KB5060954	KB5055521
KB5062799	KB5058524	KB890830

KB5054006	KB5039214	KB5029242
KB5053594	KB5037763	KB5028169
KB5052006	KB5037016	KB5027219
KB5050109	KB5036899	KB5026363
KB5049993	KB5035962	KB5025228
KB5048671	KB5035855	KB5023788
KB5046612	KB5034862	KB5023697
KB5044293	KB5034767	KB5022838
KB5043124	KB5034119	KB5022289
KB5043051	KB5033373	KB5017396
KB5041773	KB5032391	KB5012170
KB5041576	KB5032197	KB5005698
KB5040562	KB5031362	KB4589210
KB5040434	KB5030504	KB4494175
KB5039334	KB5030213	KB4132216