



Date: November 2, 2021

Author: Software Engineering

System: TRIDENT

Subsystem: Windows 7 & Software Security Patches

Subject: Windows 7 & Software Security Patch release October 2021

Purpose:

To release the list of validated Microsoft critical patches on the TRIDENT System.

Note 1:

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 001 (Windows 7 ESU License Installation) to the system prior to installation of November 2020 Windows 7 Monthly Patch Disk (SFW-03339) and April 2021 Windows 7 Monthly Patch Disk (SFW-03430) to allow for updates to be installed, due to end of support for Microsoft Windows 7.

List of patches/KB's validated for October 2021:

KB5006749	KB5006743	KB5006728
-----------	-----------	-----------

List of all patches/KB's validated on the TRIDENT System:

KB976932	KB5005615	KB5004378
KB5006749	KB5005089	KB5004307
KB5006743	KB5005088	KB5004289
KB5006728	KB5004951	KB5004120
KB5005633	KB5004757	KB5003694

KB5003667	KB4576612	KB4534314
KB5003233	KB4575903	KB4534310
KB5003228	KB4571729	KB4533012
KB5001848	KB4571719	KB4532945
KB5001392	KB4570673	KB4532932
KB5001335	KB4569775	KB4531786
KB5000851	KB4569767	KB4530734
KB5000841	KB4565623	KB4530692
KB4601363	KB4565612	KB4525235
KB4601347	KB4565539	KB4525233
KB4600945	KB4565524	KB4523206
KB4598289	KB4565354	KB4523206
KB4598279	KB4562030	KB4520003
KB4597239	KB4561669	KB4519976
KB4592510	KB4561643	KB4516655
KB4592503	KB4556843	KB4516065
KB4592471	KB4556836	KB4516033
KB4586827	KB4556798	KB4515854
KB4586805	KB4555449	KB4512506
KB4585204	KB4552940	KB4512486
KB4580970	KB4552919	KB4511516
KB4580387	KB4550965	KB4507456
KB4580345	KB4541500	KB4507449
KB4578963	KB4538483	KB4507004
KB4578952	KB4537829	KB4506997
KB4577053	KB4537813	KB4505422
KB4577051	KB4536952	KB4503292

KB4503269	KB4343899	KB3192391
KB4499175	KB4339284	KB3177467
KB4499164	KB4338823	KB3172605
KB4495606	KB4284867	KB3161949
KB4495588	KB4103712	KB3159398
KB4493472	KB4096418	KB3156016
KB4493448	KB4093108	KB3150220
KB4490628	KB4074837	KB3139914
KB4489885	KB4074587	KB3139398
KB4489878	KB4057114	KB3138910
KB4486564	KB4056897	KB3138612
KB4486563	KB4054521	KB3133977
KB4483458	KB4048960	KB3127220
KB4483451	KB4041678	KB3126587
KB4483187	KB4038779	KB3124275
KB4480970	KB4034679	KB3122648
KB4480960	KB4025337	KB3110329
KB4480063	KB4022722	KB3109560
KB4480055	KB4022619	KB3109103
KB4474419	KB4020322	KB3108664
KB4471328	KB4019990	KB3108371
KB4468323	KB4019263	KB3101722
KB4467106	KB4015546	KB3097989
KB4462915	KB4012212	KB3092601
KB4457145	KB3212642	KB3086255
KB4457044	KB3205394	KB3080079
KB4457035	KB3197867	KB3078601

KB3075220	KB2978742	KB2834140
KB3074543	KB2977292	KB2813430
KB3072305	KB2973201	KB2813347
KB3071756	KB2973112	KB2807986
KB3067903	KB2972211	KB2789645
KB3061518	KB2972100	KB2786081
KB3060716	KB2968294	KB2770660
KB3059317	KB2943357	KB2758857
KB3055642	KB2937610	KB2742599
KB3046269	KB2931356	KB2736422
KB3046017	KB2912390	KB2731771
KB3045685	KB2911501	KB2729452
KB3042058	KB2900986	KB2729094
KB3037574	KB2894854	KB2727528
KB3035126	KB2894844	KB2719033
KB3030377	KB2893294	KB2716513
KB3023215	KB2888049	KB2705219
KB3022777	KB2884256	KB2698365
KB3021674	KB2882822	KB2690533
KB3019978	KB2871997	KB2685939
KB3018238	KB2868038	KB2670838
KB3011780	KB2864202	KB2667402
KB3010788	KB2862335	KB2656356
KB2992611	KB2862330	KB2654428
KB2991963	KB2862152	KB2653956
KB2984972	KB2861698	KB2639308
KB2979596	KB2840631	KB2631813

KB2621440

KB2620704

KB2604115

KB2585542

KB2579686

KB2565063

KB2564958

KB2560656

KB2552343

KB2533623

KB2533552

KB2532531

KB2510531

KB2506212

KB2506014

KB2491683

KB2479943

KB2467173