



Date: November 9, 2022

Author: Software Engineering

System: Apex

Subsystem: Win Pro Embedded 10 LTSC 2019

Subject: Microsoft Win Pro Embedded 10 LTSC 2019 Patch release October 2022

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSC 2019 systems.

List of patches/KB's validated for October 2022:

KB5018419

List of all patches/KB's validated on the Apex System:

KB5018419	KB5012128	KB5005030
KB5017315	KB5011503	KB5003711
KB5016623	KB5010351	KB4589208
KB5015811	KB5009557	KB4583462
KB5014692	KB5008873	KB4577586
KB5014164	KB5008218	KB4535680
KB5013941	KB5007206	KB4494174
KB5013641	KB5006672	KB4022619
KB5012647	KB5005568	KB4018073
KB5012170	KB5005112	