



**Date:** November 9, 2022

**Author:** Software Engineering

**System:** TRIDENT

**Subsystem:** Windows 7 & Software Security Patches

**Subject:** Windows 7 & Software Security Patch release October 2022

---

---

**Purpose:**

To release the list of validated Microsoft critical patches on the TRIDENT System.

**Note:**

- Please check if January 2020 Windows 7 Monthly Patch Disk is installed on system, if not then first Install January 2020 Windows 7 Monthly Patch Disk (SFW-03075).
- Apply SFW-03338 Rev. 004 (Windows 7 ESU License Installation) to the system prior to installation of March 2022 Windows 7 Monthly Patch Disk (SFW-03646) to allow for updates to be installed, due to end of support for Microsoft Windows 7.

**List of patches/KB's validated for October 2022:**

KB5018479                      KB5018454

**List of all patches/KB's validated on the TRIDENT System:**

KB976932	KB5016679	KB5014742
KB5018479	KB5016676	KB5014012
KB5018454	KB5016057	KB5013999
KB5017397	KB5015862	KB5013637
KB5017373	KB5015861	KB5012649
KB5017361	KB5014748	KB5012626

KB5012138	KB5003228	KB4571729
KB5012131	KB5001848	KB4571719
KB5011649	KB5001392	KB4570673
KB5011552	KB5001335	KB4569775
KB5011529	KB5000851	KB4569767
KB5010422	KB5000841	KB4565623
KB5009621	KB4601363	KB4565612
KB5008867	KB4601347	KB4565539
KB5008859	KB4600945	KB4565524
KB5008282	KB4598289	KB4565354
KB5007233	KB4598279	KB4562030
KB5006749	KB4597239	KB4561669
KB5006743	KB4592510	KB4561643
KB5006728	KB4592503	KB4556843
KB5005633	KB4592471	KB4556836
KB5005615	KB4586827	KB4556798
KB5005089	KB4586805	KB4555449
KB5005088	KB4585204	KB4552940
KB5004951	KB4580970	KB4552919
KB5004757	KB4580387	KB4550965
KB5004378	KB4580345	KB4541500
KB5004307	KB4578963	KB4538483
KB5004289	KB4578952	KB4537829
KB5004120	KB4577053	KB4537813
KB5003694	KB4577051	KB4536952
KB5003667	KB4576612	KB4534314
KB5003233	KB4575903	KB4534310

KB4533012	KB4499164	KB4338823
KB4532945	KB4495606	KB4284867
KB4532932	KB4495588	KB4103712
KB4531786	KB4493472	KB4096418
KB4530734	KB4493448	KB4093108
KB4530692	KB4490628	KB4074837
KB4525235	KB4489885	KB4074587
KB4525233	KB4489878	KB4057114
KB4523206	KB4486564	KB4056897
KB4523206	KB4486563	KB4054521
KB4520003	KB4483458	KB4048960
KB4519976	KB4483451	KB4041678
KB4516655	KB4483187	KB4038779
KB4516065	KB4480970	KB4034679
KB4516033	KB4480960	KB4025337
KB4515854	KB4480063	KB4022722
KB4512506	KB4480055	KB4022619
KB4512486	KB4474419	KB4020322
KB4511516	KB4471328	KB4019990
KB4507456	KB4468323	KB4019263
KB4507449	KB4467106	KB4015546
KB4507004	KB4462915	KB4012212
KB4506997	KB4457145	KB3212642
KB4505422	KB4457044	KB3205394
KB4503292	KB4457035	KB3197867
KB4503269	KB4343899	KB3192391
KB4499175	KB4339284	KB3177467

KB3172605	KB3072305	KB2973201
KB3161949	KB3071756	KB2973112
KB3159398	KB3067903	KB2972211
KB3156016	KB3061518	KB2972100
KB3150220	KB3060716	KB2968294
KB3139914	KB3059317	KB2943357
KB3139398	KB3055642	KB2937610
KB3138910	KB3046269	KB2931356
KB3138612	KB3046017	KB2912390
KB3133977	KB3045685	KB2911501
KB3127220	KB3042058	KB2900986
KB3126587	KB3037574	KB2894854
KB3124275	KB3035126	KB2894844
KB3122648	KB3030377	KB2893294
KB3110329	KB3023215	KB2888049
KB3109560	KB3022777	KB2884256
KB3109103	KB3021674	KB2882822
KB3108664	KB3019978	KB2871997
KB3108371	KB3018238	KB2868038
KB3101722	KB3011780	KB2864202
KB3097989	KB3010788	KB2862335
KB3092601	KB2992611	KB2862330
KB3086255	KB2991963	KB2862152
KB3080079	KB2984972	KB2861698
KB3078601	KB2979596	KB2840631
KB3075220	KB2978742	KB2834140
KB3074543	KB2977292	KB2813430

KB2813347	KB2604115
KB2807986	KB2585542
KB2789645	KB2579686
KB2786081	KB2565063
KB2770660	KB2564958
KB2758857	KB2560656
KB2742599	KB2552343
KB2736422	KB2533623
KB2731771	KB2533552
KB2729452	KB2532531
KB2729094	KB2510531
KB2727528	KB2506212
KB2719033	KB2506014
KB2716513	KB2491683
KB2705219	KB2479943
KB2698365	KB2467173
KB2690533	
KB2685939	
KB2670838	
KB2667402	
KB2656356	
KB2654428	
KB2653956	
KB2639308	
KB2631813	
KB2621440	
KB2620704	