

Date: November 11, 2025

Author: Software Engineering

System: Apex

Subsystem: Microsoft Win Pro Embedded 10 LTSP 2016

Subject: Microsoft Win Pro Embedded 10 LTSP 2016 Patch release October 2025

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0" which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.17 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2014:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for October 2025:

KB5066836

KB5066584

List of all patches/KB's validated on the Apex System:

KB5066584	KB5039334	KB5027222
KB5066836	KB5039214	KB5027219
KB5065427	KB5037763	KB5027131
KB5065687	KB5037016	KB5026363
KB5063871	KB5036899	KB5026362
KB5062560	KB5035962	KB5025229
KB5062799	KB5035855	KB5025228
KB5061010	KB5034862	KB5023788
KB5060954	KB5034767	KB5023702
KB5058524	KB5034127	KB5023697
KB5058383	KB5034119	KB5022840
KB5055661	KB5033904	KB5022838
KB5055521	KB5033688	KB5022511
KB5054833	KB5033373	KB5022289
KB890830	KB5033371	KB5022286
KB5054006	KB5032391	KB5021237
KB5053594	KB5032197	KB5021235
KB5052006	KB5032196	KB5021126
KB5050109	KB5031984	KB5021045
KB5049993	KB5031362	KB5020866
KB5048671	KB5031361	KB5020627
KB5046612	KB5030504	KB5019966
KB5044293	KB5030214	KB5019964
KB5043124	KB5030213	KB5018419
KB5043051	KB5029931	KB5018411
KB5042749	KB5029247	KB5017396
KB5041773	KB5029242	KB5017315
KB5041576	KB5029185	KB5017305
KB5040562	KB5028169	KB5017095
KB5040434	KB5028168	KB5016884

KB5016623	KB5005573	KB4601318
KB5016622	KB5005568	KB4601060
KB5016058	KB5005112	KB4598480
KB5015811	KB5005043	KB4598243
KB5015808	KB5005030	KB4598230
KB5014702	KB5004244	KB4593226
KB5014692	KB5004238	KB4592440
KB5014553	KB5003711	KB4589210
KB5014026	KB5003646	KB4589208
KB5013952	KB5003638	KB4587735
KB5013941	KB5003243	KB4586830
KB5013641	KB5003197	KB4586793
KB5012647	KB5003171	KB4583465
KB5012596	KB5002542	KB4580346
KB5012170	KB5002520	KB4580325
KB5012128	KB5002497	KB4577668
KB5011570	KB5002464	KB4577667
KB5011503	KB5002406	KB4577586
KB5011495	KB5001407	KB4577015
KB5010359	KB5001404	KB4576750
KB5010351	KB5001402	KB4571694
KB5009557	KB5001347	KB4570333
KB5009546	KB5001342	KB4570332
KB5008873	KB5001078	KB4566424
KB5008218	KB5000859	KB4565912
KB5008207	KB5000822	KB4565511
KB5007206	KB5000803	KB4565349
KB5007192	KB5000685	KB4562562
KB5006672	KB4601393	KB4562561
KB5006669	KB4601392	KB4561616
KB5005698	KB4601345	KB4561608

KB4561600	KB4532691	KB4509091
KB4558998	KB4530715	KB4507460
KB4558997	KB4530689	KB4506998
KB4557397	KB4525236	KB4505225
KB4556813	KB4523205	KB4503308
KB4552924	KB4523204	KB4498947
KB4551853	KB4521862	KB4494440
KB4550994	KB4521858	KB4494175
KB4550929	KB4520724	KB4494174
KB4549949	KB4519998	KB4493470
KB4549947	KB4519338	KB4489882
KB4540723	KB4516115	KB4487026
KB4540670	KB4516044	KB4485447
KB4539571	KB4514366	KB4346087
KB4538461	KB4512937	KB4132216
KB4537764	KB4512578	KB4091664
KB4537759	KB4512577	KB4022619
KB4535680	KB4512574	KB4018073
KB4534273	KB4512517	KB2565063
KB4534271	KB4511553	KB2467173