



Date: November 11, 2025

Author: Software Engineering

System: BREVERA / CORLUMINA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release October

Purpose:

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.17 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for October 2025:

KB5066143

KB5066586

List of all patches/KB's validated on the BREVERA / CORLUMINA System:

KB5066143	KB5036604	KB5021126
KB5066586	KB5035849	KB5019964
KB5065428	KB5034768	KB5018411
KB5063877	KB5034127	KB5017396
KB5062557	KB5033904	KB5017305
KB5060531	KB5033371	KB5017095
KB5058392	KB5032196	KB5016884
KB5055519	KB5031984	KB5016622
KB5050533	KB5031361	KB5016058
KB890830	KB5030214	KB5015808
KB5053596	KB5029931	KB5014702
KB5052000	KB5029376	KB5014553
KB5050008	KB5029247	KB5014026
KB5049608	KB5028168	KB5013952
KB5048661	KB5027222	KB5012596
KB5046615	KB5027131	KB5012170
KB5044277	KB5026363	KB5011570
KB5044016	KB5026362	KB5011495
KB5043050	KB5025229	KB5010359
KB5042215	KB5025228	KB5009546
KB5041578	KB5023788	KB5008207
KB5040430	KB5023697	KB5007192
KB5039879	KB5022838	KB5006669
KB5039217	KB5022511	KB5005698
KB5037765	KB5022289	KB5005573
KB5036896	KB5021235	KB5005112

KB5005043	KB4571694	KB4512574
KB5004238	KB4565912	KB4512517
KB5003638	KB4565511	KB4509091
KB5003197	KB4562561	KB4507460
KB5001402	KB4561616	KB4505225
KB5001347	KB4557397	KB4498947
KB5001078	KB4556813	KB4494440
KB5000803	KB4550994	KB4494175
KB5000685	KB4550929	KB4494174
KB4601392	KB4540723	KB4493470
KB4601318	KB4540670	KB4489882
KB4598243	KB4537764	KB4487026
KB4593226	KB4535680	KB4485447
KB4589210	KB4534271	KB4346087
KB4589208	KB4530689	KB4132216
KB4586830	KB4525236	KB4091664
KB4580346	KB4521858	KB2565063
KB4577586	KB4520724	KB2467173
KB4577015	KB4519998	
KB4576750	KB4516044	