



**Date:** November 11, 2025

**Author:** Software Engineering

**System:** Path Vision

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release October 2025

---

---

**Purpose:**

To release the list of validated Microsoft critical patches on the Path Vision System.

Note: Core Vision and BioVision are on similar product line as Path Vision and both use the same executable so patches validated on Path Vision system can be installed on Core Visions system and BioVision as well.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.17 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for October 2025:**

KB5066586                      KB5066143

**List of all patches/KB's validated on the Path Vision System:**

KB5066143	KB5062557	KB890830
KB5066586	KB5060531	KB5053596
KB5065428	KB5058392	KB5052000
KB5063877	KB5055519	KB5050008

KB5049608	KB5033371	KB5020627
KB5048661	KB5032196	KB5019966
KB5046615	KB5031984	KB5018419
KB5044277	KB5031361	KB5017315
KB5044016	KB5030214	KB5016623
KB5043050	KB5029931	KB5015811
KB5041578	KB5029247	KB5014692
KB5040430	KB5028168	KB5013941
KB5039879	KB5027222	KB5013641
KB5039217	KB5027131	KB5012170
KB5037765	KB5026362	KB5005112
KB5036896	KB5025229	KB4589208
KB5036604	KB5023702	KB4577586
KB5035849	KB5022840	KB4535680
KB5034768	KB5022511	KB4494174
KB5034127	KB5021237	
KB5033904	KB5020866	