



Date: November 11, 2025

Author: Software Engineering

System: SELENIA DIMENSIONS / 3DIMENSIONS

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release October 2025

Purpose:

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.17 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for October 2025:

KB5066836

KB5066584

**List of all patches/KB's validated on the SELENIA DIMENSIONS /
3DIMENSIONS System:**

KB5066584	KB5066836	KB5065427
KB5065687	KB5040562	KB5023697
KB5063871	KB5040434	KB5022838
KB5062560	KB5039334	KB5022289
KB5062799	KB5039214	KB5021235
KB5060954	KB5037763	KB5021126
KB5061010	KB5037016	KB5019964
KB5058524	KB5036899	KB5018411
KB5058383	KB5035962	KB5017396
KB5055661	KB5035855	KB5017305
KB5055521	KB5034862	KB5017095
KB5050533	KB5034767	KB5016884
KB890830	KB5034119	KB5016622
KB5054006	KB5033373	KB5016058
KB5053594	KB5032391	KB5015808
KB5052006	KB5032197	KB5014702
KB5050109	KB5031362	KB5014553
KB5049993	KB5030504	KB5014026
KB5048671	KB5030213	KB5013952
KB5046612	KB5029376	KB5013625
KB5044293	KB5029242	KB5012596
KB5043124	KB5028169	KB5012170
KB5043051	KB5027219	KB5012118
KB5042215	KB5026363	KB5011570
KB5041773	KB5025228	KB5011495
KB5041576	KB5023788	KB5010359

KB5009546	KB4576750	KB4503267
KB5008877	KB4571694	KB4498947
KB5008207	KB4565912	KB4494440
KB5007192	KB4565511	KB4494175
KB5006669	KB4562561	KB4493472
KB5005698	KB4561616	KB4493470
KB5005573	KB4557397	KB4493448
KB5005043	KB4556813	KB4493435
KB5004238	KB4550994	KB4489882
KB5003638	KB4550929	KB4485447
KB5003197	KB4540723	KB4346087
KB5001402	KB4540670	KB4132216
KB5001347	KB4537764	KB4091664
KB5001078	KB4535680	KB2565063
KB5000803	KB4534271	KB2467173
KB5000685	KB4530689	
KB4601392	KB4525236	
KB4601318	KB4521858	
KB4601051	KB4520724	
KB4598243	KB4519998	
KB4593226	KB4516044	
KB4589210	KB4512574	
KB4586830	KB4512517	
KB4580346	KB4509091	
KB4578969	KB4507460	
KB4577586	KB4505225	
KB4577015	KB4503537	