



Date: November 11, 2025

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2016 (v1607) & Software Security Patches

Subject: Windows 10 LTSC 2016 (v1607) & Software Security Patch release October 2025

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.17 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for October 2025:

KB5066836

KB5066584

List of all patches/KB's validated on the SecurView System:

KB5066584

KB5062799

KB5055521

KB5066836

KB5061010

KB890830

KB5065427

KB5060954

KB5054006

KB5065687

KB5058524

KB5053594

KB5063871

KB5058383

KB5052006

KB5062560

KB5055661

KB5050109

KB5049993	KB5028169	KB5008207
KB5048671	KB5027219	KB5007192
KB5046612	KB5026363	KB5006669
KB5044293	KB5025228	KB5005698
KB5043124	KB5023788	KB5005573
KB5043051	KB5023697	KB5005112
KB5041773	KB5022838	KB5005043
KB5041576	KB5022286	KB5004238
KB5040562	KB5021235	KB5003638
KB5040434	KB5020866	KB5003197
KB5039334	KB5019964	KB5001402
KB5039214	KB5018411	KB5001347
KB5037763	KB5017396	KB5001078
KB5037016	KB5017305	KB5000803
KB5036899	KB5017095	KB4601393
KB5035962	KB5016622	KB4601392
KB5035855	KB5016058	KB4601345
KB5034862	KB5015808	KB4601318
KB5034767	KB5014702	KB4601060
KB5034119	KB5014026	KB4598243
KB5033373	KB5013952	KB4593226
KB5032391	KB5012596	KB4589210
KB5032197	KB5012170	KB4589208
KB5031362	KB5011570	KB4586830
KB5030504	KB5011495	KB4580346
KB5030213	KB5010359	KB4580325
KB5029242	KB5009546	KB4577586

KB4577015	KB4512574
KB4576750	KB4512517
KB4571694	KB4509091
KB4565912	KB4503308
KB4565511	KB4494175
KB4562561	KB4494174
KB4561616	KB4346087
KB4561600	KB4132216
KB4557397	KB2565063
KB4556813	KB2467173
KB4550994	
KB4550949	
KB4550929	
KB4540723	
KB4540670	
KB4539571	
KB4537764	
KB4537759	
KB4535680	
KB4534271	
KB4530689	
KB4525236	
KB4521858	
KB4520724	
KB4519998	
KB4516115	
KB4516044	