



Date: November 11, 2025

Author: Software Engineering

System: Vision DXA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release October 2025

Purpose:

To release the list of validated Microsoft critical patches on the Vision DXA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.17 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for October 2025:

KB5066143

KB5066586

List of all patches/KB's validated on the Vision DXA System:

KB5066143

KB5055519

KB5046615

KB5066586

KB890830

KB5044016

KB5065428

KB5053596

KB5041578

KB5063877

KB5052000

KB5040430

KB5062557

KB5050008

KB5039879

KB5060531

KB5049608

KB5039217

KB5058392

KB5048661

KB5037765

KB5036896	KB5029247	KB5020627
KB5036604	KB5028168	KB5019966
KB5035849	KB5027222	KB5018419
KB5034768	KB5027131	KB5017315
KB5034127	KB5026362	KB5016623
KB5033904	KB5025229	KB5015811
KB5033371	KB5023702	KB5013641
KB5032196	KB5022840	KB5012170
KB5031984	KB5022511	KB5005112
KB5031361	KB5022286	KB4589208
KB5030214	KB5021237	KB4577586
KB5029931	KB5020866	KB4535680