



Date: October 21, 2024

Author: Software Engineering

System: SecurView

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
September 2024

Purpose:

To release the list of validated Microsoft critical patches on the SecurView System.

List of patches/KB's validated for September 2024:

KB5043050

List of all patches/KB's validated on the SecurView System:

KB890830	KB5034119	KB5027222
KB5043050	KB5033904	KB5027131
KB5041578	KB5033371	KB5026362
KB5040430	KB5032391	KB5025229
KB5039879	KB5032196	KB5023702
KB5039217	KB5031984	KB5022840
KB5037765	KB5031361	KB5022511
KB5036896	KB5030214	KB5022286
KB5036604	KB5029931	KB5021237
KB5035849	KB5029376	KB5021126
KB5034768	KB5029247	KB5020866
KB5034127	KB5028168	KB5020627

KB5019966	KB5012647	KB5007206
KB5018419	KB5012170	KB5005698
KB5017315	KB5012128	KB5005112
KB5016623	KB5011503	KB4589210
KB5015811	KB5010351	KB4589208
KB5014692	KB5009557	KB4577586
KB5013941	KB5008873	KB4535680
KB5013641	KB5008218	KB4132216