



Date: October 10, 2025

Author: Software Engineering

System: Cenova

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
September 2025

Purpose:

To release the list of validated Microsoft critical patches on the Cenova System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for September 2025:

KB5065428

List of all patches/KB's validated on the Cenova System:

KB5065428	KB5060531	KB5050533
KB5063877	KB5058392	KB890830
KB5062557	KB5055519	KB5053596

KB5052000	KB5029247	KB5012128
KB5050008	KB5028168	KB5011503
KB5049608	KB5027222	KB5010351
KB5048661	KB5027131	KB5009557
KB5046615	KB5027124	KB5008873
KB5044277	KB5026362	KB5008218
KB5044016	KB5025229	KB5007206
KB5043050	KB5023702	KB5006672
KB5042215	KB5022840	KB5005568
KB5041578	KB5022511	KB5005112
KB5040430	KB5022286	KB5005030
KB5039879	KB5021237	KB5004244
KB5039217	KB5021126	KB5003711
KB5037765	KB5020866	KB5003646
KB5036896	KB5020627	KB5003638
KB5036604	KB5019966	KB5003243
KB5035849	KB5018419	KB5003171
KB5034768	KB5017315	KB5001407
KB5034127	KB5016884	KB5001404
KB5033904	KB5016623	KB5001402
KB5033371	KB5015811	KB5001342
KB5032196	KB5014692	KB5000859
KB5031984	KB5014553	KB5000822
KB5031361	KB5013941	KB5000685
KB5030214	KB5013641	KB4601393
KB5029931	KB5012647	KB4601345
KB5029376	KB5012170	KB4601060

KB4589210

KB4577586

KB4535680

KB4589208

KB4557397

KB4494174

KB4580325

KB4539571

KB4465065