



**Date:** October 10, 2025

**Author:** Software Engineering

**System:** SELENIA DIMENSIONS / 3DIMENSIONS

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release  
September 2025

---

**Purpose:**

To release the list of validated Microsoft critical patches on the SELENIA DIMENSIONS / 3DIMENSIONS System.

**Latest .NET Versions Validated for the Product:**

**Note:** Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

**Latest SQLSERVER Versions Validated:**

**Note:** Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

**List of patches/KB's validated for September 2025:**

KB5065428

**List of all patches/KB's validated on the SELENIA DIMENSIONS /  
3DIMENSIONS System:**

KB5065428	KB5036604	KB5020615
KB5063877	KB5035849	KB5019966
KB5062557	KB5034768	KB5018419
KB5060531	KB5034127	KB5017315
KB5058392	KB5033904	KB5016884
KB5055519	KB5033371	KB5016623
KB5050533	KB5032196	KB5015811
KB5049615	KB5031990	KB5014692
KB890830	KB5031361	KB5014553
KB5053596	KB5030214	KB5014164
KB5052000	KB5029925	KB5013941
KB5050008	KB5029376	KB5013641
KB5049608	KB5029247	KB5013626
KB5048661	KB5028168	KB5012647
KB5046615	KB5027222	KB5012170
KB5044277	KB5027124	KB5012119
KB5044016	KB5026362	KB5011503
KB5043050	KB5025229	KB5010351
KB5042215	KB5023702	KB5009557
KB5041578	KB5022840	KB5008873
KB5040430	KB5022504	KB5008218
KB5039879	KB5022286	KB5007206
KB5039217	KB5021237	KB5006672
KB5037765	KB5021126	KB5005568
KB5036896	KB5020874	KB5005112

KB5005030	KB4566424
KB5004244	KB4565349
KB5003711	KB4562562
KB5003646	KB4561608
KB5003243	KB4561600
KB5003171	KB4558998
KB5001407	KB4558997
KB5001404	KB4557397
KB5001342	KB4552924
KB5000859	KB4551853
KB5000822	KB4549949
KB5000685	KB4549947
KB4601393	KB4539571
KB4601345	KB4537759
KB4601060	KB4535680
KB4598480	KB4514366
KB4598230	KB4505225
KB4592440	KB4494174
KB4589208	KB4022619
KB4587735	KB4018073
KB4586793	KB2565063
KB4580325	KB2467173
KB4577668	
KB4577667	
KB4577586	
KB4570333	
KB4570332	