



Date: October 10, 2025

Author: Software Engineering

System: SecurView Manager – Windows Server 2016

Subsystem: Windows 10 Server 2016 & Software Security Patches

Subject: Windows 10 Server 2016 & Software Security Patch release September 2025

Purpose:

To release the list of validated Microsoft critical patches on the SecurView Manager – Windows Server 2016 System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for September 2025:

KB5065427

KB5065687

List of all patches/KB's validated on the SecurView Manager – Windows Server 2016 System:

KB5065427

KB5065687

KB5063871

KB5058524

KB5053594

KB5062560

KB5055661

KB5052006

KB5062799

KB5055521

KB5050109

KB5061010

KB890830

KB5049993

KB5060954

KB5054006

KB5048671

KB5046612	KB5035962	KB5027219
KB5044293	KB5035855	KB5026363
KB5043124	KB5034862	KB5025228
KB5043051	KB5034767	KB5023788
KB5041773	KB5034119	KB5023697
KB5041576	KB5033373	KB5022838
KB5040562	KB5032391	KB5022289
KB5040434	KB5032197	KB5017396
KB5039334	KB5031362	KB5012170
KB5039214	KB5030504	KB5005698
KB5037763	KB5030213	KB4589210
KB5037016	KB5029242	KB4494175
KB5036899	KB5028169	KB4132216