



Date: October 10, 2025

Author: Software Engineering

System: Vision DXA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release
September 2025

Purpose:

To release the list of validated Microsoft critical patches on the Vision DXA System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.14 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for September 2025:

KB5065428

List of all patches/KB's validated on the Vision DXA System:

KB5065428	KB5053596	KB5041578
KB5063877	KB5052000	KB5040430
KB5062557	KB5050008	KB5039879
KB5060531	KB5049608	KB5039217
KB5058392	KB5048661	KB5037765
KB5055519	KB5046615	KB5036896
KB890830	KB5044016	KB5036604

KB5035849	KB5027222	KB5018419
KB5034768	KB5027131	KB5017315
KB5034127	KB5026362	KB5016623
KB5033904	KB5025229	KB5015811
KB5033371	KB5023702	KB5013641
KB5032196	KB5022840	KB5012170
KB5031984	KB5022511	KB5005112
KB5031361	KB5022286	KB4589208
KB5030214	KB5021237	KB4577586
KB5029931	KB5020866	KB4535680
KB5029247	KB5020627	
KB5028168	KB5019966	