

**HOLOGIC®**

**HTx40 Scope**  
Compatible with 5Fr  
instrumentation

**HTx60s Scope**  
Compatible with MyoSure®  
Lite, Reach and Manual devices

**HTx2000  
Image Processing System**



# HTx

## Single-Use Scope System

Introducing the newest innovation in the Hologic hysteroscopic portfolio

The HTx System features two single-use scope options for a wide range of hysteroscopic procedures.

# See and treat with the new HTx Single-Use Scope System

## EASE OF USE


Portable system that allows you to easily capture, manage, and export case footage.

## PERFORMANCE

Precise, quality image with autofocus and white balance.

## VALUE

All-in-one single-use system designed for the office and operating room.



### HTx40 vs. HT60s comparison

Specifications	HTx40	HTx60s
Outer diameter (OD)	4.5 mm	6.2 mm
Working channel	2.1 mm	3.3 mm
Working length	240 mm	195 mm
Direction of view	8°	12°

### Ordering information

Product description	Product number	QTY
<b>HYSTEROSCOPES</b>		
HTx40 Scope (4.5 OD)	HT21003	6/box
HTx60s Scope (6.2 OD)	HT21007	6/box
<b>CAPITAL</b>		
HTx2000 Image Processing System (System includes one image cable)	HT21072	1
<b>INSTRUMENTS</b>		
MyoSure Manual Device	20-403ML	3/box
MyoSure Lite Device	30-403LITE	3/box
Single Action Scissors, Reusable	HT21010	1
Biopsy Forceps, Reusable	HT21028	1
Alligator Grasper, Reusable	HT21046	1
Hook Grasper, Reusable	HT21064	1

Product description	Product number	QTY
<b>ACCESSORIES</b>		
Procedure Pack	HT22014	6/box
<b>OPTIONAL</b>		
Omni 4K UHD Office Monitor	83-12-4001	1
Image Cable Kit	HT22007	3/box

Note: The image processing unit serves as a portable monitor, so a secondary monitor may not be needed.

#### References:

HTx Disposable Hysteroscope System IFU

### Contact your Hologic sales representative for more information.

All Rights Reserved. Hologic, MyoSure, Omni and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. HTx is manufactured by Suzhou Acuvu Medical Technology Co., Ltd. and distributed exclusively in the United States by Hologic.

Hologic, Inc. | 250 Campus Dr., Marlborough, MA 01752 | 508 263 2900  
 SS-01488-001 Rev 002 | 04/2024 | © 2024 Hologic, Inc.



Hologic.com