

Date: August 07, 2026

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2016 & Software Security Patches

Subject: Windows Server 2016 & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5099542

KB5099535

KB5090354

List of all patches/KB's validated on the AWM System:

KB5099542

KB5002879

KB5094141

KB5099535

KB5087065

KB5087537

KB5090354

KB5094122

KB5088064

KB5082089	KB5049993	KB5028169
KB5082198	KB5048671	KB5027219
KB5068402	KB5046612	KB5026363
KB5078938	KB5044293	KB5025228
KB5075999	KB5043124	KB5023788
KB5075902	KB5043051	KB5023697
KB5073447	KB5042215	KB5022838
KB5073722	KB5041773	KB5022289
KB5071543	KB5041578	KB5021235
KB5065225	KB5041576	KB5021126
KB5068864	KB5040562	KB5019964
KB5070247	KB5040434	KB5018411
KB5066584	KB5039879	KB5017396
KB5066836	KB5039334	KB5017305
KB5065687	KB5039214	KB5017095
KB5065427	KB5037763	KB5016884
KB5063871	KB5037016	KB5016622
KB5062560	KB5036899	KB5016058
KB5062799	KB5035962	KB5015808
KB5060954	KB5035855	KB5014702
KB5061010	KB5034862	KB5014553
KB5058392	KB5034767	KB5014026
KB5049608	KB5034119	KB5013952
KB5055661	KB5033373	KB5012596
KB5055521	KB5032391	KB5012170
KB5050533	KB5032197	KB5011570
KB890830	KB5031362	KB5011495
KB5054006	KB5030504	KB5010359
KB5053594	KB5030213	KB5009557
KB5052006	KB5029376	KB5009546
KB5050109	KB5029242	KB5008873

KB5008207	KB4557397
KB5007192	KB4551853
KB5006669	KB4550929
KB5005698	KB4549947
KB5005573	KB4540723
KB5005112	KB4539571
KB5005043	KB4537759
KB5004238	KB4535680
KB5003638	KB4514366
KB5003197	KB4505225
KB5001402	KB4494175
KB5001347	KB4494174
KB5001078	KB4132216
KB5000803	KB4132214
KB5000685	KB2565063
KB4601392	KB2467173
KB4601318	
KB4598243	
KB4593226	
KB4589210	
KB4589208	
KB4586830	
KB4580346	
KB4577586	
KB4577015	
KB4576750	
KB4571694	
KB4565912	
KB4565511	
KB4562561	
KB4561616	