

Date: September 8, 2026

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2019 & Software Security Patches

Subject: Windows Server 2019 & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

<u>Date</u>	<u>List of KB's</u>		
July 2026	KB5100989	KB5090354	KB5099538
August 2026	KB5120698	KB5120238	

List of all patches/KB's validated on the AWM System before July 2026:

KB5100989	KB5042215	KB5019964
KB5090354	KB5041578	KB5018419
KB5099538	KB5040430	KB5018411
KB5094123	KB5039879	KB5017396
KB5087061	KB5039217	KB5017305
KB5087538	KB5037765	KB5017095
KB5068402	KB5036896	KB5016884
KB5082123	KB5036604	KB5016622
KB5082413	KB5035849	KB5016058
KB5078752	KB5034768	KB5015808
KB5075904	KB5034127	KB5014702
KB5073723	KB5033904	KB5014553
KB5071544	KB5033371	KB5014026
KB5065225	KB5032196	KB5013952
KB5068791	KB5031984	KB5013641
KB5066143	KB5031361	KB5012596
KB5066586	KB5030214	KB5012170
KB5049608	KB5029931	KB5011570
KB5063877	KB5029376	KB5011495
KB5062557	KB5029247	KB5010359
KB5060531	KB5028168	KB5009557
KB5058392	KB5027222	KB5009546
KB5055519	KB5027131	KB5008873
KB5050533	KB5026362	KB5008207
KB890830	KB5025229	KB5007192
KB5053596	KB5023702	KB5006669
KB5052000	KB5022840	KB5005698
KB5050008	KB5022511	KB5005573
KB5049608	KB5022286	KB5005112
KB5048661	KB5021237	KB5005043
KB5046615	KB5021126	KB5004238
KB5044277	KB5020866	KB5003638
KB5044016	KB5020627	KB5003197
KB5043050	KB5019966	KB5001402
KB5001347	KB4589210	KB4565912
KB5001078	KB4589208	KB4565511
KB5000803	KB4586830	KB4562561
KB5000685	KB4580346	KB4561616
KB4601392	KB4577586	KB4557397
KB4601318	KB4577015	KB4551853
KB4598243	KB4576750	KB4550929
KB4593226	KB4571694	KB4549947
KB4540723	KB4505225	KB2467173
KB4539571	KB4494175	KB2565063
KB4537759	KB4494174	KB4514366
KB4535680	KB4132216	