

Date: September 08, 2026

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2022 & Software Security Patches

Subject: Windows Server 2022 & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's:

<u>Date</u>	<u>KB's validated:</u>		
July 2026	KB5101010	KB5099540	KB5090354
August 2026	KB5120705	KB5120242	

List of all patches/KB's validated on the AWM System before July 2026:

KB5094128	KB5087068	KB5068402
-----------	-----------	-----------