

Date: August 07, 2026

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2022 & Software Security Patches

Subject: Windows Server 2022 & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5101010
KB5099540
KB5090354

List of all patches/KB's validated on the AWM System:

KB5101010
KB5099540
KB5090354

KB5094128

KB5087068

KB5068402