

Date: July 06, 2026

Author: Software Engineering

System: AWM

Subsystem: Windows Server 2022 & Software Security Patches

Subject: Windows Server 2022 & Software Security Patch release June 2026

Purpose:

To release the list of validated Microsoft critical patches on the AWM System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

- **.NET Core 2.1 / 2.2**

Hosting Bundle 2.1.30

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for June 2026:

KB5094128

KB5087068

KB5068402

List of all patches/KB's validated on the AWM System:

KB5094128

KB5087068

KB5068402