

Date: September 04, 2026

Author: Software Engineering

System: Affirm Prone

Subsystem: Windows 10 & Software Security Patches

Subject: Windows 10 & Software Security Patch release August 2026

Purpose:

To release the list of validated Microsoft critical patches on the Affirm Prone System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated:

<u>Date</u>	<u>List of KB's</u>		
July 2026	KB5100989	KB5099538	KB5090354
August 2026	KB5120698	KB5120238	

List of all patches/KB's validated on the AFFIRM PRONE System before July 2026:

KB5100989	KB5068402	KB5065428
KB5099538	KB5078752	KB5063877
KB5090354	KB5075904	KB5062557
KB5094123	KB5073723	KB5060531
KB5087061	KB5065225	KB5058392
KB5087538	KB5068791	KB5055519
KB5082413	KB5066143	KB5050533
KB5082123	KB5066586	KB890830
KB5053596	KB5027131	KB5005568
KB5052000	KB5026362	KB5005112
KB5050008	KB5025229	KB5005030
KB5049608	KB5023702	KB5004244
KB5048661	KB5022840	KB5003711
KB5046615	KB5022511	KB5003646
KB5044277	KB5022286	KB5003243
KB5044016	KB5021237	KB5003171
KB5043050	KB5021126	KB5001407
KB5042215	KB5020866	KB5001404
KB5041578	KB5020627	KB5001342
KB5040430	KB5019966	KB5000859
KB5039879	KB5018419	KB5000822
KB5039217	KB5017315	KB5000685
KB5037765	KB5016884	KB4601393
KB5036896	KB5016623	KB4601345
KB5036604	KB5015811	KB4601060
KB5035849	KB5014692	KB4598480
KB4598230	KB4558998	KB4521862
KB4592440	KB4558997	KB4519338
KB4589208	KB4557397	KB4516115
KB4587735	KB4552924	KB4514366
KB4586793	KB4551853	KB4512937
KB4580325	KB4549949	KB4512578
KB4577668	KB4549947	KB4512577
KB4577667	KB4539571	KB4511553
KB4577586	KB4538461	KB4506998
KB4570333	KB4537759	KB4505225
KB4570332	KB4535680	KB4503308

KB4566424	KB4534273	KB4494174
KB4565349	KB4532691	KB2565063
KB4562562	KB4530715	KB2467173
KB4561608	KB4523205	KB4561600
KB4523204		