

Date: August 07, 2026

Author: Software Engineering

System: Affirm Prone

Subsystem: Windows 10 & Software Security Patches

Subject: Windows 10 & Software Security Patch release July 2026

Purpose:

To release the list of validated Microsoft critical patches on the Affirm Prone System.

Latest .NET Versions Validated for the Product:

Note: Security updates should only be installed for existing .NET installations already present on the product. For example, .NET 8.0.24 update can be installed only if an earlier version of .NET 8.0 is already part of the product installation.

- **.NET Framework 3.5 / 4.x**

The latest security updates are included in the Windows Cumulative Update.

Latest SQLSERVER Versions Validated:

Note: Security updates should only be installed for existing SQLSERVER installations already present on the product.

- **SQL Server 2017:**

The latest security updates are included in the Windows Cumulative Update.

List of patches/KB's validated for July 2026:

KB5100989	KB5099538	KB5090354
-----------	-----------	-----------

List of all patches/KB's validated on the AFFIRM PRONE System:

KB5100989	KB5068402	KB5065428
KB5099538	KB5078752	KB5063877
KB5090354	KB5075904	KB5062557
KB5094123	KB5073723	KB5060531
KB5087061	KB5065225	KB5058392
KB5087538	KB5068791	KB5055519
KB5082413	KB5066143	KB5050533
KB5082123	KB5066586	KB890830

KB5053596	KB5027131	KB5005568
KB5052000	KB5026362	KB5005112
KB5050008	KB5025229	KB5005030
KB5049608	KB5023702	KB5004244
KB5048661	KB5022840	KB5003711
KB5046615	KB5022511	KB5003646
KB5044277	KB5022286	KB5003243
KB5044016	KB5021237	KB5003171
KB5043050	KB5021126	KB5001407
KB5042215	KB5020866	KB5001404
KB5041578	KB5020627	KB5001342
KB5040430	KB5019966	KB5000859
KB5039879	KB5018419	KB5000822
KB5039217	KB5017315	KB5000685
KB5037765	KB5016884	KB4601393
KB5036896	KB5016623	KB4601345
KB5036604	KB5015811	KB4601060
KB5035849	KB5014692	KB4598480
KB5034768	KB5014553	KB4598230
KB5034127	KB5013941	KB4592440
KB5033904	KB5013641	KB4589208
KB5033371	KB5012647	KB4587735
KB5032196	KB5012170	KB4586793
KB5031984	KB5012128	KB4580325
KB5031361	KB5011503	KB4577668
KB5030214	KB5010351	KB4577667
KB5029931	KB5009557	KB4577586
KB5029376	KB5008873	KB4570333
KB5029247	KB5008218	KB4570332
KB5028168	KB5007206	KB4566424
KB5027222	KB5006672	KB4565349

KB4562562	KB4538461	KB4514366
KB4561608	KB4537759	KB4512937
KB4561600	KB4535680	KB4512578
KB4558998	KB4534273	KB4512577
KB4558997	KB4532691	KB4511553
KB4557397	KB4530715	KB4506998
KB4552924	KB4523205	KB4505225
KB4551853	KB4523204	KB4503308
KB4549949	KB4521862	KB4494174
KB4549947	KB4519338	KB2565063
KB4539571	KB4516115	KB2467173